



# ThermaLite® Disposable Electrodes for the TCU®!

## Disposable Convenience *plus* a Faster Heat-Up Time!

A nice addition to the TCU accessories! Fits handles for both Geiger TCU and original Geiger Cautery Model 100. Conveniently packaged in boxes of 10 for easy access and storage. Non-sterile.



- Quality, disposable palladium tips
- Hollow brass shanks for faster heat-up
- Available in three styles: 210T Loop, 212 Angled Loop and 214 Fine tip

Order by contacting your Geiger Medical Technologies Distributor or Geiger Medical Technologies.



### TCU THERMALITE® PLATINUM DISPOSABLE ELECTRODES

Item # 210T - Loop Tip  
212 - Angled Loop Tip  
214 - Fine Tip



[www.geigermedical.com](http://www.geigermedical.com)

**GEIGER**  
medical technologies

608 13th Avenue, Council Bluffs, IA 51501 USA  
Toll Free: 1.888.GEIGER.1 (1.888.434.471) or 712.323.3432 • Fax: 800.320.9612 or 712.323.1156 • Email: [info@geigermedical.com](mailto:info@geigermedical.com)  
Geiger is a registered trademark of Geiger Medical Technologies. ©2003 Geiger Medical Technologies. #LT 03-0017